

IHE Integration Statement

Orion Health is committed to an open standards approach to information exchange. All of our products are designed to use technologies and standards endorsed and in use by healthcare organizations, thereby allowing integration with systems built by other vendors as part of a wider integrated clinical information system.

IHE promotes the coordinated use of established standards to address specific needs in the support of optimal patient care.

VENDOR	SYSTEM NAME	SYSTEM VERSION	DATE
Orion Health	Health Information Exchange	2	January 1, 2013

INTEGRATION PROFILE IMPLEMENTED	ACTOR IMPLEMENTED	INTEGRATION PROFILE OPTION IMPLEMENTED
[ATNA] Audit Trail and Node Authentication	Secure Application	None
[CT] Consistent Time	Time Client	None
[MPQ] Multi-Patient Queries	Document Consumer	None
[PDQ] Patient Demographics Query	Patient Demographics Consumer	Patient Demographics and Visit Query
[PDQ] Patient Demographics Query	Patient Demographics Consumer	On-Demand Documents
[PDQ] Patient Demographics Query	Patient Demographics Supplier	None
[PDQ] Patient Demographics Query	PDS-Patient Demographics Supplier	None
[PDQv3] Patient Demographics Query HL7 V3	Patient Demographics Supplier	None
[PIX] Patient Identifier Cross-referencing for MPI	Patient Identifier Cross-reference Consumer	PIX Update Notification
[PIX] Patient Identifier Cross-referencing for MPI	Patient Identifier Cross-reference Consumer	None
[PIX] Patient Identifier Cross-referencing for MPI	Patient Identifier Cross-reference Manager	None
[PIXv3] Patient Identifier Cross-reference HL7 V3	Patient Identifier Cross-reference Manager	None
[XCA] Cross-Community Access	Initiating Gateway	None
[XCA] Cross-Community Access	Responding Gateway	None
[XCPD] Cross-Community Patient Discovery	Initiating Gateway	None

[XCPD] Cross-Community Patient Discovery	Responding Gateway	None
[XDR] Cross-Enterprise Document Reliable Interchange	Document Recipient	None
[XDR] Cross-Enterprise Document Reliable Interchange	Document Recipient	Accepts Limited Metadata
[XDR] Cross-Enterprise Document Reliable Interchange	Document Source	None
[XDS] Cross-Enterprise Document Sharing	Document Consumer	None
[XDS] Cross-Enterprise Document Sharing	Document Consumer	On-Demand Documents
[XDS] Cross-Enterprise Document Sharing	Document Registry	None
[XDS] Cross-Enterprise Document Sharing	Document Registry	On-Demand Documents
[XDS] Cross-Enterprise Document Sharing	Document Respository	None
[XDS] Cross-Enterprise Document Sharing	Document Source	Folder Management
[XDS] Cross-Enterprise Document Sharing	Document Source	None
[XDS] Cross-Enterprise Document Sharing	Document Source	Document Addendum
[XDS] Cross-Enterprise Document Sharing	Document Source	Document Transformation
[XDS] Cross-Enterprise Document Sharing	Document Source	Document Replace
[XDS] Cross-Enterprise Document Sharing	On-Demand Document Source	Persistence of Retrieved Documents
[XDS] Cross-Enterprise Document Sharing	On-Demand Document Source	None
[XPHR] Exchange of Personal Health Record Content	Content Creator	None

INTERNET ADDRESS FOR VENDOR'S IHE INFORMATION

<https://orionhealth.com/us/knowledge-hub/disclosures/>

LINKS TO STANDARDS CONFORMANCE STATEMENTS FOR THE IMPLEMENTATION HL7 DICOM

LINKS TO GENERAL INFORMATION ON THE IHE

In North America www.ihe.net
 In Europe www.ihe-europe.org
 In Japan www.jira-net.or.jp/ihe-j

Find out more at orionhealth.com