

IHE Integration Statement

VENDOR	SYSTEM NAME	SYSTEM VERSION	DATE
Orion Health	Rhapsody	5	January 4, 2012

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

INTEGRATION PROFILE IMPLEMENTED	ACTOR IMPLEMENTED	INTEGRATION PROFILE OPTION IMPLEMENTED
[ATNA] Audit Trail and Node Authentication	Secure Application	None
[CT] Consistent Time	Time Client	None
[MPQ] Multi-Patient Queries	Document Consumer	None
[PDQ] Patient Demographics Query	Patient Demographics Consumer	Patient Demographics and Visit Query
[PDQ] Patient Demographics Query	Patient Demographics Consumer	None
[PIX] Patient Identifier Cross-referencing for MPI	Patient Identifier Cross-reference Consumer	PIX Update Notification
[PIX] Patient Identifier Cross-referencing for MPI	Patient Identifier Cross-reference Consumer	None
[XDR] Cross-Enterprise Document Reliable Interchange	Document Recipient	None
[XDR] Cross-Enterprise Document Reliable Interchange	Document Recipient	Accepts Limited Metadata
[XDR] Cross-Enterprise Document Reliable Interchange	Document Source	None
[XDS] Cross-Enterprise Document Sharing	Document Consumer	None
[XDS] Cross-Enterprise Document Sharing	Document Consumer	On-Demand Documents
[XDS] Cross-Enterprise Document Sharing	Document Repository	None
[XDS] Cross-Enterprise Document Sharing	Document Source	Folder Management
[XDS] Cross-Enterprise Document Sharing	Document Source	None

[XDS] Cross-Enterprise Document Sharing	Document Source	Document Addendum
[XDS] Cross-Enterprise Document Sharing	Document Source	Document Transformation
[XDS] Cross-Enterprise Document Sharing	Document Source	Document Replace

INTERNET ADDRESS FOR VENDOR'S IHE INFORMATION

<https://orionhealth.com/us/knowledge-hub/disclosures/>

LINKS TO STANDARDS CONFORMANCE STATEMENTS FOR THE IMPLEMENTATION

HL7

DICOM

LINKS TO GENERAL INFORMATION ON THE IHE

In North America	www.ihe.net
In Europe	www.ihe-europe.org
In Japan	www.jira-net.or.jp/ihe-j

Find out more at orionhealth.com

Rhapsody® Integration Engine is intended only for the electronic transfer, storage, or display of medical device data, or the electronic conversion of such data from one format to another in accordance with a preset specification as specified in the product manual and/or related documentation. Rhapsody Integration Engine is not intended to be used for active patient monitoring, controlling or altering the functions or parameters of any medical device, or any other purpose relating to data obtained directly or indirectly from a medical device other than the transfer, storage, and conversion of such data from one format to another in accordance with preset specifications. Orion Health makes no warranties and the functionality described within may change without notice. Rhapsody® is a registered trademark of Orion Health™ Limited, manufactured in New Zealand, by Orion Health Limited. All other trademarks displayed in this document are the property of Orion Health or their respective owners, and may not be used without written permission of the owner. Rhapsody Integration Engine is not intended to be used for diagnostic purposes, or to replace clinical judgment or responsibilities.